



ProbeHunter Can Now Test All OEMs' Ultrasound Probes

Posted by Melanie Hamilton-Basich | Feb 3, 2021 | Testing & Calibration | 0



BBS Medical, Stockholm, Sweden, has released an adapter that allows its ProbeHunter to conduct dynamic testing and validation of the performance of Mindray DC8 series ultrasound probes, just the latest addition to the company's probe testing capabilities.

"This means that we now can test all leading OEMs' ultrasound probes—GE, Philips, Siemens, Samsung, Mindray, and more—independent from the system, on ProbeHunter," says BBS Medical CEO Björn Segall. "With 50 different multi-brand adapters in our portfolio, we can offer hospitals and anyone dealing or servicing ultrasound probes the ultimate tool for patient safety."

ProbeHunter is a multi-brand real-time testing device with innovative technology designed to provide quality assurance for ultrasound probes. It is designed for patient safety, for the user to get a super quick "Go or No Go" before using a probe on a patient.

"It is now easier than ever to get a ProbeHunter system," adds Segall. "The Swedish government supports export business and is offering a capital credit guarantee during the pandemic. It means, in general, that it is possible to have a ProbeHunter system today and pay later. I would like to encourage hospitals and the ultrasound probe industry to make use of this unique possibility for patient safety."

Featured image: CEO Björn Segall testing a GE-4c-RS ultrasound probe using ProbeHunter. (Courtesy: BBC Medical/ProbeHunter)

SHARE:



RECENT POSTS

[Pylo Launches 4G Blood Pressure Monitor for RPM](#)

[ProbeHunter Can Now Test All OEMs' Ultrasound Probes](#)

[AdaptHealth Acquires Home Medical Equipment Provider AeroCare](#)

[Inovytect's Ventway Sparrow Ventilators Nab FDA Clearance](#)

[The Right to Repair Movement Is Poised to Explode in 2021](#)

FOLLOW US

[Facebook](#)

[Twitter](#)

[LinkedIn](#)

DELIVER THE RIGHT MESSAGE, TO THE RIGHT PERSON, ON THE RIGHT CHANNEL, ON THE RIGHT DEVICE, AT THE RIGHT TIME.

*The complex sale is dead.
Only the complex purchase exists.*

KAIROS
INTELLIGENCE-BASED
MARKETING

[< PREVIOUS](#)

[AdaptHealth Acquires Home Medical Equipment Provider AeroCare](#)

[NEXT >](#)

[Pylo Launches 4G Blood Pressure Monitor for RPM](#)

RELATED POSTS





Medical Test and Measurement Equipment Market to See Growth

November 27, 2017



Testing Device for Electrosurgical Generators

April 15, 2013



Rigel Electrical Analyzers Selected for South Africa Project

March 28, 2012



Fluke Biomedical Offers New Infusion Device Analyzer

April 22, 2014

Leave a reply

Your email address will not be published. Required fields are marked *

COMMENT

NAME *

EMAIL *

WEBSITE

☐ Save my name, email, and website in this browser for the next time I comment.

POST COMMENT

SEARCH ...

FOLLOW US



Facebook



Twitter



LinkedIn

MEDQOR MEDIA BRANDS

24x7 Magazine

Axis Imaging News

Clinical Lab Products

Orthodontic Products

Physical Therapy Products

Plastic Surgery Practice

Rehab Management

RT Magazine

Sleep Review

The Hearing Review

RESOURCES

24x7 Magazine

Videos

Webinars

White Papers

HELPFUL LINKS

Media Kit

Subscribe

Contact Us