

2023-06-15

## ProbeHunter Patent approved Europe!



We are happy to announce that ProbeHunter technology is now patented in Europe as well as in USA and Sweden.

With more than 65 different multi-brand adapters in our portfolio, we can offer hospitals and anyone dealing or servicing ultrasound probes the ultimate tool for patient safety and QA.

ProbeHunter is a multi-brand real-time dynamic testing device with innovative technology designed to provide quality assurance for ultrasound probes. It is designed for patient safety, for the user to get a super quick “Go or No Go” before using a probe on a patient.

BBS Medical, Stockholm, Sweden, has released new adapters that allows its ProbeHunter to conduct dynamic testing and validation of the performance of all leading Endoscopic Ultrasound that utilizes advanced features for US-imaging, sampling and interventional EUS/EBUS procedures this includes Pentax, Olympus and FUJI EUS probes. This is the latest addition to the company’s probe testing capabilities.

This means that we can include all leading EUS OEM’s to the adapter portfolio and to the adapter list of all OEMs’ ultrasound probes—GE, Philips, Canon, Siemens, Samsung, Mindray, and more—to test independent from the US-system, on ProbeHunter.

Feel free to contact me for more info and how we can support your hospital or business!

BR

Bjorn

[www.ProbeHunter.com](http://www.ProbeHunter.com)