TechNation

YOU ARE AT: Home » Magazine » Industry Updates » News and Notes » Europe Approves Probehunter Patent

Europe Approves Probehunter Patent

JUNE 19, 2023 🕚 1 M N READ

ProbeHunter technology is now patented in Europe as well as in USA and Sweden.

With more than 65 different multi-brand adapters in its portfolio, ProbeHunter can offer hospitals and anyone dealing or servicing ultrasound probes the "ultimate tool for patient safety and QA," according to a release.

ProbeHunter is a multi-brand real-time dynamic testing device with innovative technology designed to provide quality assurance for ultrasound probes. It is designed for patient safety, for the user to get a quick "Go or No Go" before using a probe on a patient. BBS Medical in Stockholm, Sweden, has released new adapters that allows its ProbeHunter to conduct dynamic testing and validation of the performance of all leading endoscopic ultrasound that utilizes advanced features for ultrasound-imaging, sampling and interventional EUS/EBUS procedures this includes Pentax, Olympus and FUJI EUS probes. This is the latest addition to the company's probe testing capabilities.

This means that ProbeHunter can include all leading EUS OEMs to the adapter portfolio and to the adapter list of all OEMs' ultrasound probes — GE, Philips, Canon, Siemens, Samsung, Mindray and more — to test independent from the ultrasound system, on ProbeHunter.

For more information, visit ProbeHunter.com.

f

 \square